

## Product Datasheet

### Z-VAD-FMK (orb1147018)

<b>Description</b>	Pan-caspase inhibitor; Peptides.
<b>Form/Appearance</b>	White lyophilised solid
<b>Storage</b>	Store dry, dark and frozen
<b>Note</b>	For research use only
<b>Protein Sequence</b>	Z-VAD-FMK
<b>Formula</b>	$C_{22}H_{30}FN_3O_7$
<b>Purity</b>	> 98% by TLC
<b>MW</b>	467.49 Da
<b>CAS Number</b>	187389-52-2
<b>Solubility (25°C)</b>	Soluble to 10 mg/mL in DMSO
<b>Expiration Date</b>	12 months from date of receipt.