

Product Datasheet

TRAP (orb1147027)

Description	Thrombin receptor agonist; Peptides.
Form/Appearance	Freeze dried solid
Storage	Store dry, frozen and in the dark
Note	For research use only
Protein Sequence	SFLLRNPNDKYEPF
Formula	$C_{81}H_{118}N_{20}O_{23}$
Purity	> 95% by hplc
MW	1739.94 Da
CAS Number	137339-65-2
Solubility (25°C)	Soluble to 0.90 mg/ml in water
Expiration Date	12 months from date of receipt.