

Product Datasheet

SFFLRN-NH2 (orb1147042)

Description	Protease activated receptor-1 (PAR-1) activator; Peptides.
Form/Appearance	Freeze dried solid
Storage	Store desiccated, frozen and in the dark
Note	For research use only
Protein Sequence	SFFLRN-NH2
Formula	$C_{37}H_{55}N_{11}O_8$
Purity	> 95% by HPLC
MW	781.9 Da
Solubility (25°C)	Soluble in water to 1mg per ml
Expiration Date	12 months from date of receipt.