

Product Datasheet

Oxyntomodulin (orb1147065)

Description	Endogenous glucagon like peptide receptor agonist; Peptides.
Form/Appearance	Freeze dried solid
Storage	Store dry, dark and frozen
Note	For research use only
Protein Sequence	HSQGTFTSDYSKYLDSRRAQDFVQWLMNTKRNKNNIA
Formula	$C_{192}H_{295}N_{59}O_{60}S$
Purity	> 95% by HPLC
MW	4421.86 Da
CAS Number	62340-29-8
Solubility (25°C)	Soluble in dilute acid
Expiration Date	12 months from date of receipt.