

## Product Datasheet

### Camostat mesylate (orb1147116)

<b>Description</b>	Serine protease inhibitor; Small molecules.
<b>Form/Appearance</b>	Off-white powder
<b>Storage</b>	Store desiccated, frozen and in the dark
<b>Note</b>	For research use only
<b>Protein Sequence</b>	4-2-2-dimethylamino-2-oxoethoxy-2-oxoethylphenyl 4-guanidinobenzoate, methanesulfonic acid salt
<b>Formula</b>	$C_{20}H_{22}N_4O_5 \cdot CH_3SO_3H$
<b>Purity</b>	> 98% by HPLC
<b>MW</b>	494.5 Da
<b>CAS Number</b>	59721-29-8
<b>Solubility (25°C)</b>	Soluble in water and DMSO to 50mg/ml
<b>Expiration Date</b>	12 months from date of receipt.