

## Product Datasheet

### Rabbit IgG control (orb1147195)

<b>Description</b>	Polyclonal antibody to Rabbit IgG Isotype control
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC, LF
<b>Preservatives</b>	PBS
<b>Form/Appearance</b>	Liquid/ Lyophilized
<b>Storage</b>	Store in -20°C
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Protein A purified
<b>Expiration Date</b>	12 months from date of receipt.