

## Product Datasheet

### E2 CSFV Monoclonal antibody (Swine) (orb1147213)

<b>Description</b>	E2 CSFV Monoclonal antibody (Swine)
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	Recombinant Swine CSFV E2 Protein
<b>Preservatives</b>	1X PBS + 0.001% rAlbumin pH 7.9
<b>Form/Appearance</b>	Liquid/ Lyophilized
<b>Storage</b>	Store in -20°C
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Protein A purified
<b>MW</b>	150kda
<b>Expiration Date</b>	12 months from date of receipt.