

Product Datasheet

Zebrafish Caveolin 3/CAV3 Antibody (orb1147633)

Description	Zebrafish Caveolin 3/CAV3 Antibody
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Unconjugated
Form/Appearance	Liquid
Concentration	500µg/ml
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	Q6YLN7
Expiration Date	12 months from date of receipt.