



Product Datasheet

JAM-A/F11R Antibody (orb1147774)


Description	JAM-A/F11R Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IF, IHC, WB
Immunogen	E.coli-derived human JAM-A/F11R recombinant protein (Position: S28-T294).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	38 kDa
Uniprot ID	Q9Y624
Expiration Date	12 months from date of receipt.



IHC analysis of JAM-A/F11R using anti-JAM-A/F11R antibody



IHC analysis of JAM-A/F11R using anti-JAM-A/F11R antibody



IF analysis of JAM-A/F11R using anti-JAM-A/F11R antibody