



Product Datasheet

CD74 Antibody (orb1147813)

Description	CD74 (the MHC II-associated invariant chain, Ii) is a type II transmembrane protein expressed in ant...
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, IP, WB
Immunogen	SU-DHL-4 lymphoma cells
Target	CD74
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 2-5 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P04233
Entrez	972
Dilution Range	Flow cytometry: Recommended dilution: 2-5 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human CD74 positive lympho...



Flow cytometry surface staining pattern ...