

---

## Product Datasheet

### ECFP-Tag Monoclonal Antibody (6B11) (orb1147816)

Evolution Biocomponents.  
Description

Featured Anti-ECFP-Tag Monoclonal

	your tagged protein quantif...
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Recombinant Protein
<b>Preservatives</b>	Liquid in PBS, pH 7.4, containing 0.02% Sodium Azide as preservative and 50% Glycerol.
<b>Form/Appearance</b>	Liquid solution
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:2000-1:5000).
<b>Isotype</b>	Mouse IgG
<b>Clonality</b>	Monoclonal
<b>Hazard Information</b>	The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.
<b>Expiration Date</b>	12 months from date of receipt.