
Product Datasheet

MBP Tag Mouse Monoclonal Antibody (9Y5) (orb1147837)

Evolution Bioreagents.
Description

Featured Anti-MBP Tag Mouse

	designed for your tagged protein qua...
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Preservatives	Liquid in PBS, pH 7.4, containing 0.02% Sodium Azide as preservative and 50% Glycerol.
Form/Appearance	Liquid solution
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:5000).
Isotype	Mouse IgG
Clonality	Monoclonal
Hazard Information	The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.
Expiration Date	12 months from date of receipt.