

---

## Product Datasheet

### Agarose Conjugated Anti-Myc Tag Mouse Monoclonal Antibody (2D5) (orb1147840)

**Description**

## Featured Anti-Myc Tag Mouse

specially designed for immunoprecipitation...

**Species/Host**

Mouse

**Conjugation**

Agarose

**Tested Applications**

IP

**Immunogen**

Synthetic Peptide

**Preservatives**

50% gel slurry suspended in PBS, pH 7.4, containing 0.02% Sodium Azide as preservative.

**Form/Appearance**

Liquid solution

**Storage**

Store at 4°C for 12 months. Avoid freeze-thaw or centrifugation.

**Note**

For research use only

**Application notes**

Optimal working dilutions should be determined experimentally by the investigator.

**Isotype**

Mouse IgG1

**Clonality**

Monoclonal

**Hazard Information**

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.

**Expiration Date**

12 months from date of receipt.