
Product Datasheet

HRP Conjugated Anti- β -Tubulin Mouse Monoclonal Antibody (3G6) (orb1147860)

DescriptionFeatured Anti- β -Tubulin Mouse

Conjugated, specially designed for you...

Species/Host

Mouse

Reactivity

Canine, Gallus, Hamster, Human, Insect, Monkey, Mouse, Rabbit, Rat, Sheep, Yeast

Conjugation

HRP

Tested Applications

IHC-P, WB

Immunogen

Synthetic Peptide

Preservatives

Liquid in PBS, pH 7.4 and 50% Glycerol.

Form/Appearance

Liquid solution

Storage

Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

Note

For research use only

Application notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:5000), IHC-P (1:200).

Isotype

Mouse IgG1

Clonality

Monoclonal

Hazard Information

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.