

## Product Datasheet

### Anti-Valproic Acid Antibody (orb1148009)

<b>Description</b>	Anti-Valproic Acid Mouse Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Valproic acid-KLH.
<b>Target</b>	Valproic Acid
<b>Preservatives</b>	PBS, pH 7.4.
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified
<b>CAS Number</b>	99-66-1
<b>Expiration Date</b>	12 months from date of receipt.