

Product Datasheet

Anti-Valproic Acid Mouse Monoclonal Antibody - Library Pack (orb1148010)

Description	Anti-Valproic Acid Mouse Monoclonal Antibody Library Pack - Includes 100 ug of each of the following...
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Valproic acid conjugated to KLH.
Target	Valproic Acid
Preservatives	PBS, pH 7.4.
Form/Appearance	Liquid
Storage	This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Purity	Purified
CAS Number	99-66-1
Expiration Date	12 months from date of receipt.