

Product Datasheet

Anti-Transmissible Gastroenteritis Virus Antibody (orb1148023)

Description	Anti-Transmissible Gastroenteritis Mouse Monoclonal Antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, NeA
Immunogen	TGEV-infected ST cells
Target	Transmissible Gastroenteritis Virus
Preservatives	PBS, pH 7.4
Form/Appearance	Liquid
Storage	This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Dilution Range	These antibodies have been qualified for use in ELISA to detect TGEV. #18003 and #18004 may be used in TGEV neutralization assays. Endusers should determine optimal concentrations for their applications.
Expiration Date	12 months from date of receipt.