

## Product Datasheet

# Anti-Transmissible Gastroenteritis Virus Antibody (orb1148025)

<b>Description</b>	Anti-Transmissible Gastroenteritis Mouse Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	TGEV-infected ST cells
<b>Target</b>	Transmissible Gastroenteritis Virus
<b>Preservatives</b>	PBS, pH 7.4
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified
<b>Dilution Range</b>	These antibodies have been qualified for use in ELISA to detect TGEV. #18003 and #18004 may be used in TGEV neutralization assays. Endusers should determine optimal concentrations for their applications.
<b>Expiration Date</b>	12 months from date of receipt.