
Product Datasheet

Anti-Transmissible Gastroenteritis Mouse

Monoclonal Antibody - Library Pack (orb1148027)

Description

Anti-Transmissible

Monoclonal Antibody Library
Pack - Includes 100 ug of each ...**Species/Host**

Mouse

Reactivity

Virus

Conjugation

Unconjugated

**Tested
Applications**

ELISA

ImmunogenTGEV-infected ST cells Host
Species: Mouse Specificity:
These antibodies recognize
TGEV and do not cross-react
with FIPV, FeCV, or CCV. #18003
and #18004 can neutralize TGEV
in vitro.**Target**

Transmissible Gastroenteritis

Preservatives

PBS, pH 7.4

Form/Appearance

Liquid

StorageStore frozen product at or below
-20C. Thawed product may be
stored for 2-4 weeks at 4C. For
optimal storage, aliquot and
store at -20C. Avoid multiple
freeze- thaw cycles.**Note**

For research use only

Clonality

Monoclonal

Purity

Purified

Dilution RangeThese antibodies have been
qualified for use in ELISA to
detect TGEV. #18003 and
#18004 may be used in TGEV
neutralization assays. Endusers
should determine optimal
concentrations for their
applications.**Expiration Date**

12 months from date of receipt.