

## Product Datasheet

### Anti-Theophylline Antibody (orb1148031)

<b>Description</b>	Anti-Theophylline Mouse Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Theophylline-KLH.
<b>Target</b>	Theophylline
<b>Preservatives</b>	PBS, pH 7.4.
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified
<b>CAS Number</b>	58-55-9
<b>Dilution Range</b>	These antibodies have been qualified for use in ELISA to detect theophylline.
<b>Expiration Date</b>	12 months from date of receipt.