

Product Datasheet

Anti-TAG-72 Antibody (orb1148048)

Description	Anti-TAG-72 Mouse Monoclonal Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant human TAG-72.
Target	TAG-72
Preservatives	0.1% Sodium Azide. Tris buffer, pH 7.3-7.7, 1% rAlbumin, 0.1% sodium azide.
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at 2-8C.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Dilution Range	Immunohistochemistry: use at a dilution of 1:100-1:200 on formalin-fixed, paraffin-embedded samples after heat-induced epitope retrieval at pH 9 for 10-30 minutes.
Expiration Date	12 months from date of receipt.