

## Product Datasheet

### Anti-SARS Matrix Protein (orb1148095)

<b>Description</b>	Anti-SARS Matrix Protein Rabbit Polyclonal Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Synthetic peptide corresponding to the amino-terminus of SARS-CoV matrix protein.
<b>Target</b>	SARS Matrix Protein
<b>Preservatives</b>	0.02% Sodium Azide. PBS, 0.02% sodium azide.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1.0 mg/mL
<b>Storage</b>	This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Purified
<b>Expiration Date</b>	12 months from date of receipt.