

Product Datasheet

Anti-SARS Envelope Protein (orb1148096)

Description	Anti-SARS Envelope Protein Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Synthetic peptide corresponding to the carboxy-terminus of SARS-CoV envelope protein.
Target	SARS Envelope Protein
Preservatives	0.02% Sodium Azide. PBS, 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1.0 mg/mL
Storage	This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.