

Product Datasheet

Anti-SARS Envelope Protein (orb1148097)

| | |
|----------------------------|---|
| Description | Anti-SARS Envelope Protein Rabbit Polyclonal Antibody |
| Species/Host | Rabbit |
| Reactivity | Virus |
| Conjugation | Unconjugated |
| Tested Applications | ELISA |
| Immunogen | Synthetic peptide corresponding to the amino-terminus of SARS-CoV envelope protein. |
| Target | SARS Envelope Protein |
| Preservatives | 0.02% Sodium Azide. PBS, 0.02% sodium azide. |
| Form/Appearance | Liquid |
| Concentration | 1.0 mg/mL |
| Storage | This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| Purity | Purified |
| Expiration Date | 12 months from date of receipt. |