

Product Datasheet

Anti-Rabbit Immunoglobulin Antibody (orb1148140)

Description	Anti-Rabbit Immunoglobulin Mouse Monoclonal Antibody
Species/Host	Mouse
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Pool of purified rabbit immunoglobulin from normal sera
Target	Rabbit Immunoglobulin
Preservatives	PBS, pH 7.4.
Form/Appearance	Liquid
Storage	These antibodies are stable for at least one (1) year at -20C to -70C. Store product in appropriate aliquots to avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Dilution Range	These antibodies have been qualified for use in ELISA to detect rabbit immunoglobulin.
Expiration Date	12 months from date of receipt.