
Product Datasheet

Anti-Rabbit Immunoglobulin Mouse Monoclonal Antibody - Library Pack (orb1148141)

Description

Anti-Rabbit Immunoglobulin

Library Pack - Includes 100 ug of each of the f...

Species/Host

Mouse

Reactivity

Rabbit

Conjugation

Unconjugated

Tested

ELISA

Applications
Immunogen

Pool of purified rabbit immunoglobulin from normal sera. Host Species: Mouse
 Specificity: #14601 recognizes rabbit IgG and IgM; #14602 recognizes rabbit IgG. These antibodies do not cross-react with immunoglobulin of the following species: bovine, canine, chicken, donkey, equine, feline, rat, sheep, porcine, guinea pig, goat, or human.

Target

Rabbit Immunoglobulin

Preservatives

PBS, pH 7.4.

Form/Appearance

Liquid

Storage

Store at 2-8C. Shelf life is one (1) year from date of receipt.

Note

For research use only

Clonality

Monoclonal

Purity

Purified

Dilution Range

These antibodies have been qualified for use in ELISA to detect rabbit immunoglobulin.

Expiration Date

12 months from date of receipt.