
Product Datasheet

**Anti-Rabbit IgG H&L Chain Antibody, Cross
Adsorbed, HRP conjugate (orb1148144)**

Description

Anti-Rabbit IgG, H&L Chain

	Polyclonal Antibody
Species/Host	Goat
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC, WB
Immunogen	Rabbit IgG Host Species: Goat Stabilizers: 0.2% Bovine Serum Albumin (BSA) Preservatives: 0.1% Pro-Clean 400 Buffer: Phosphate-buffered saline (PBS) Purification: Immunoaffinity chromatography
Target	Rabbit IgG H&L Chain, Cross Adsorbed, HRP conjugate
Preservatives	0.1% Pro-Clean 400. 0.2% Bovine Serum Albumin (rAlbumin), 0.1% Pro-Clean 400, Phosphate-buffered saline (PBS)
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store at 2-8C. Shelf life is one (1) year from date of receipt.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Dilution Range	ELISA: 1:10,000-1:50,000 Western Blot: 1:500-1:30,000 Immunohistochemistry: 1:100-1:250 Immunocytochemistry: 1:100-1:250 Prepare working dilutions immediately before use. These are recommended dilutions. User should determine optimal dilutions for their applications.
Expiration Date	12 months from date of receipt.