
Product Datasheet

Anti-Rabbit IgG Fc specific Antibody, FITC conjugate (orb1148148)

Description

Anti-Rabbit IgG, Fc Specific-FITC Goat

Species/Host

Goat

Reactivity

Rabbit

Conjugation

FITC

Tested

ICC, IHC

Applications**Immunogen**

Rabbit IgG Fc fragment Host Species:
Goat Stabilizers: 0.2% Bovine Serum
Albumin (BSA) Preservatives: 0.09%
Sodium Azide Buffer: Phosphate-
buffered saline (PBS) Purification:
Immunoaffinity chromatography
Conjugation: F/P ratio is 3.7

Target

Rabbit IgG Fc specific, FITC conjugate

Preservatives

0.09% Sodium Azide. 0.2% Bovine
Serum Albumin (rAlbumin), 0.09%
Sodium Azide, Phosphate-buffered
saline (PBS)

Form/Appearance

Liquid

Concentration

1.0 mg/mL

Storage

Store at 2-8C. Shelf life is one (1)
year from date of receipt.

Note

For research use only

Clonality

Polyclonal

Purity

Purified

Dilution Range

Immunohistochemistry: 1:200-
1:2,000 Immunocytochemistry:
1:200-1:2,000 Prepare working
dilutions immediately before use.
These are recommended dilutions.
User should determine optimal
dilutions for their applications.

Expiration Date

12 months from date of receipt.