

Product Datasheet

Anti-Norovirus/Calicivirus Antibody (orb1148190)

Description	Anti-Norovirus/Calicivirus Mouse Monoclonal Antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Rabbit hemorrhagic disease virus + feline calicivirus
Target	Norovirus/Calicivirus
Preservatives	0.09% Sodium Azide. PBS, pH 7.4, 0.09% sodium azide.
Form/Appearance	Liquid
Storage	This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Dilution Range	Sandwich ELISA: Recommended pairs: Human norovirus detection 18928 capture - 18929 detect 18929 capture - 18931 detect Feline calicivirus detection 18927 capture - 18931 detect 18928 capture - 18931 detect 18929 capture - 18929 detect Endusers should determine optimal antibody concentrations relative to their virus preparations.
Expiration Date	12 months from date of receipt.