

## Product Datasheet

### Anti-Norovirus/Calicivirus Antibody (orb1148192)

<b>Description</b>	Anti-Norovirus/Calicivirus Mouse Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Rabbit hemorrhagic disease virus + feline calicivirus
<b>Target</b>	Norovirus/Calicivirus
<b>Preservatives</b>	0.09% Sodium Azide. PBS, pH 7.4, 0.09% sodium azide.
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	This antibody is stable for at least one (1) year at -20C. Avoid repeated freezing and thawing.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified
<b>Dilution Range</b>	Sandwich ELISA: Recommended pairs: Human norovirus detection 18928 capture - 18929 detect 18929 capture - 18931 detect Feline calicivirus detection 18927 capture - 18931 detect 18928 capture - 18931 detect 18929 capture - 18929 detect
<b>Expiration Date</b>	12 months from date of receipt.