
Product Datasheet

Anti-NFkB-p65 (Phospho-Ser468) Antibody (orb1148197)

Description

Anti-NFkB-p65 (Phospho-Ser468)

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Peptide sequence that includes phosphorylation site of serine 468 derived from human NF κ B-p65 and conjugated to KLH. Accession no.: Q04206, NP_001138610.1 Host Species: Rabbit Specificity: This antibody detects endogenous human, mouse, and rat NF κ B-p65 only when phosphorylated at serine 468.
Target	NFkB-p65 (Phospho-Ser468)
Preservatives	0.02% Sodium Azide. PBS (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1.0 mg/mL
Storage	This antibody is stable for at least one (1) year at -20C. Can be stored at 4C for short-term use.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Dilution Range	Western blotting: use at dilution of 1:500- 1:1,000. A band of ~65kDa is detected. Immunofluorescence: use at dilution of 1:100- 1:200. These are recommended working dilutions. Enduser should determine optimal dilutions for their applications.
Expiration Date	12 months from date of receipt.