
Product Datasheet

Anti-NFkB-p65 (Phospho-Ser276) Antibody (orb1148198)

Description

Anti-NFkB-p65 (Phospho-Ser276) Rabbit

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Peptide sequence that includes phosphorylation site of serine 276 (R-P-S(p)-D-R) derived from human NF κ B-p65 and conjugated to KLH. Accession no.: Q04206, NP_001138610.1 Host Species: Rabbit Specificity: This antibody detects endogenous human, mouse, and rat NF κ B-p65 only when phosphorylated at serine 276.
Target	NFkB-p65 (Phospho-Ser276)
Preservatives	0.02% Sodium Azide. PBS (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1.0 mg/mL
Storage	This antibody is stable for at least one (1) year at -20C. Can be stored at 4C for short-term use.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Dilution Range	Western blotting: use at dilution of 1:500-1:1,000. A band of ~65kDa is detected. Immunohistochemistry: use at dilution of 1:50-1:100. These are recommended working dilutions. Enduser should determine optimal dilutions for their applications.
Expiration Date	12 months from date of receipt.