
Product Datasheet

Anti-Mouse IgG Fc specific Antibody, Biotin conjugate (orb1148278)

Description

Anti-Mouse IgG Fc Fragment-Biotin

Species/Host	Goat
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA, ICC, IHC, WB
Immunogen	<p>Mouse IgG-Fc fragment Host Species: Goat Specificity: By immunoblotting and ELISA this antibody reacts specifically with mouse IgG. Cross-reactivity with mouse IgA and IgM is negligible. This antibody may cross-react with IgG of other species.</p>
Target	Mouse IgG Fc specific, Biotin conjugate
Preservatives	0.09% Sodium Azide. PBS (PBS), 0.2% rAlbumin, 0.09% sodium azide, conjugated to biotin.
Form/Appearance	Liquid
Concentration	1.0 mg/mL
Storage	This antibody is stable for at least 1 year at 2-8C.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Dilution Range	<p>Immunoblotting: use at a dilution of 1:10,000-1:200,000.</p> <p>Immunohistochemistry: use at a dilution of 1:250-1:2,500.</p> <p>Immunocytochemistry: use at a dilution of 1:100-1:500. ELISA: use at a dilution of 1:10,000-1:200,000. These are recommended dilutions. Endusers should determine optimal dilutions for their applications. Prepare working dilutions immediately before use.</p>
Expiration Date	12 months from date of receipt.