
Product Datasheet

Anti-Mouse IgG Fc specific Antibody (orb1148280)

Description

Anti-Goat Mouse IgG, Fc Goat Polyclonal

Species/Host	Goat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC, WB
Immunogen	Mouse IgG Fc fragment Host Species: Goat Stabilizers: None Preservatives: 0.1% Sodium Azide Buffer: Phosphate-buffered saline (PBS), pH 7.2 Purification: Immunoaffinity chromatography Conjugation: None
Target	Mouse IgG Fc specific
Preservatives	0.1% Sodium Azide. 0.1% Sodium Azide, Phosphate-buffered saline (PBS), pH 7.2
Form/Appearance	Liquid
Concentration	1.0 mg/mL
Storage	This antibody is stable for at least 1 year at 2-8C.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Dilution Range	ELISA: 1:1,000-1:30,000 Western Blot: 1:1,000-1:30,000 Immunohistochemistry: 1:200-1:2,000 Immunocytochemistry: 1:200-1:2,000 Prepare working dilutions immediately before use. These are recommended dilutions. User should determine optimal dilutions for their applications.
Expiration Date	12 months from date of receipt.