

---

## Product Datasheet

### Anti-Human Metapneumovirus (hMPV) Antibody (orb1148332)

**Description**

Anti-Human Metapneumovirus (hMPV)

<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IF, WB
<b>Immunogen</b>	hMPV. Host species: Mouse. Antibody Class: IgG1. Preservative: None.
<b>Target</b>	Human Metapneumovirus (hMPV)
<b>Preservatives</b>	Lyophilized, 0.1M Tris, 0.1M glycine, 2% sucrose, 1mg/ml.
<b>Form/Appearance</b>	Lyophilized
<b>Concentration</b>	1.0 mg/mL
<b>Storage</b>	This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified
<b>Dilution Range</b>	ELISA: use at 0.1-1.0ug/ml. Immunoblotting: use at 1-5ug/ml. Immunofluorescence: use at 2-10ug/ml on virus-infected cells. These are recommended concentrations. Enduser should determine optimal concentrations for their applications.
<b>Expiration Date</b>	12 months from date of receipt.