

Product Datasheet

Anti-HIV-1 gp24 Antibody (orb1148361)

Description	Anti-HIV-1 gp24 Mouse Monoclonal Antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	r.gp24 Host Species: Mouse Subtype: IgG1 Clone Number: HIV1gp24 Stabilizers: None Preservatives: None
Target	HIV-1 gp24
Preservatives	PBS, pH 7.4, Lyophilized
Form/Appearance	Lyophilized
Storage	This antibody is stable for at least one year. Store lyophilized product at 4C. After reconstitution, aliquot and store at -20C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.