

## Product Datasheet

### Anti-Hepatitis E Virus Antibody (orb1148407)

<b>Description</b>	Anti-Hepatitis E Virus Rabbit Polyclonal Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Peptide corresponding to amino acids 627-637 of the 71 kD Hepatitis E Virus (HEV) capsid protein. Host Species: Rabbit Stabilizers: None Preservatives: 0.02% sodium azide.
<b>Target</b>	Hepatitis E Virus
<b>Preservatives</b>	PBS, pH 7.4
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	This antibody is stable for at least one (1) year at -70C. Avoid multiple freeze- thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Purified
<b>Dilution Range</b>	Immunoblotting: use at 1:500-1:1,000.
<b>Expiration Date</b>	12 months from date of receipt.