

Product Datasheet

Anti-Hepatitis C Virus Antibody (orb1148409)

Description	Anti-Hepatitis C Virus Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide corresponding to amino acids 9-21 of Hepatitis C Virus (HCV) core protein C1. Host Species: Rabbit Stabilizers: None Preservatives: 0.02% sodium azide.
Target	Hepatitis C Virus
Preservatives	PBS, pH 7.4
Form/Appearance	Liquid
Storage	These antibodies are stable for at least one (1) year at -20C to -70C. Store product in appropriate aliquots to avoid multiple freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Dilution Range	Immunoblotting: use at 1:500-1:1,000.
Expiration Date	12 months from date of receipt.