

---

## Product Datasheet

### Anti-Goat IgG H&L Chain Antibody, HRP conjugate (orb1148431)

**Description**

Anti-Goat IgG H&amp;L Chain, Peroxidase

<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Goat
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	ELISA, ICC, IHC, WB
<b>Immunogen</b>	Goat IgG Host Species: Rabbit Stabilizers: None Preservatives: None Buffer: Phosphate-buffered saline (PBS) containing 0.1% Pro-Clean 400 Purification: Immunoaffinity chromatography Conjugation: Molar enzyme:antibody ratio is 4:1
<b>Target</b>	Goat IgG H&L Chain, HRP conjugate
<b>Preservatives</b>	0.1% Pro-Clean 400. Phosphate-buffered saline (PBS) containing 0.1% Pro-Clean 400
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1.0 mg/mL
<b>Storage</b>	Store at 2-8C. Shelf life is one (1) year from date of receipt.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Purified
<b>Dilution Range</b>	ELISA: 1:10,000-1:100,000 Western Blot: 1:1,000-1:30,000 Immunohistochemistry: 1:200-1:500 Immunocytochemistry: 1:200-1:500 Prepare working dilutions immediately before use. These are recommended dilutions. User should determine optimal dilutions for their applications.
<b>Expiration Date</b>	12 months from date of receipt.