

Product Datasheet

Anti-FK-506 Antibody (orb1148445)

Description	Anti-FK-506 Mouse Monoclonal Antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	FK-506 conjugated to bovine serum albumin. Host Species: Mouse Mouse Isotype: IgM Stabilizers: None Preservatives: 0.1% sodium azide.
Target	FK-506
Preservatives	0.1% Sodium Azide. PBS containing 0.1% sodium azide
Form/Appearance	Liquid
Storage	This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
Purity	Purified
Dilution Range	This antibody may be used in ELISA and Western blotting to detect FK-506.
Expiration Date	12 months from date of receipt.