
Product Datasheet

Anti-Otoferlin Antibody (orb1148448)

Description

Anti-Otoferlin Mouse Monoclonal Antibody

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	GST fusion protein corresponding to aa 1-395 of human otoferlin (accession no. Q9HC10). Host species: Mouse. Antibody Class: IgG1. Preservative: None.
Target	Otoferlin
Preservatives	Lyophilized, 0.1M Tris, 0.1M glycine, 2% sucrose, 1mg/ml.
Form/Appearance	Lyophilized
Concentration	1.0 mg/mL
Storage	This product is stable for at least one (1) year at -20C to -70C. Reconstituted product should be stored in appropriate aliquots to avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	Q9HC10
Dilution Range	Immunoblotting: use at 1-5ug/ml. A band of ~220kDa is detected. Immunofluorescence: use at 2-10ug/ml. These are recommended concentrations. Enduser should determine optimal concentrations for their applications.
Expiration Date	12 months from date of receipt.