

Product Datasheet

Anti-PID/MTA2 Antibody (orb1148658)

Description	Anti-PID/MTA2 Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Peptide corresponding to aa 652-668 of human PID/MTA2. This sequence differs from that of mouse by one amino acid. Host Species: Rabbit Stabilizers: None Preservatives: 0.02% sodium azide
Target	PID/MTA2
Preservatives	PBS, pH 7.4.
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -20C. Avoid multiple freeze- thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	O94776
Dilution Range	Immunoblotting: use at 0.5-1 ug/ml. In immunoblots, a band of 75 kD is detected. Positive control: HeLa cell lysate.
Expiration Date	12 months from date of receipt.