

Product Datasheet

Anti-Canine parvovirus Antibody (orb1148897)

Description	Anti-Canine parvovirus Mouse Monoclonal Antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Canine parvovirus
Target	Canine parvovirus
Preservatives	0.09% Sodium Azide. PBS, pH 7.4, 0.09% sodium azide.
Form/Appearance	Liquid
Storage	This antibody is stable for 1 month at 4°C or one year at -20°C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Dilution Range	TABLE~~~~~
Expiration Date	12 months from date of receipt.