

---

## Product Datasheet

### Anti-Influenza B Neuraminidase Antibody (orb1148909)

**Description**

Anti-Influenza B Neuraminidase Mouse

<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, NeA
<b>Immunogen</b>	19601 = Influenza B Lee/40; 19602, 19603 = Influenza B HK/73. Host Species: Mouse Antibody Class: IgG2a Stabilizers: None. Preservatives: None. Available on request.
<b>Target</b>	Influenza B Neuraminidase
<b>Preservatives</b>	PBS, pH 7.4.
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	These antibodies are stable for at least one (1) year at -20 to -70°C. Store product in appropriate aliquots to avoid multiple freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified
<b>Uniprot ID</b>	<b>P03474</b>
<b>Dilution Range</b>	These antibodies may be used in immunoassays to detect and quantitate Influenza B neuraminidase and in neuraminidase inhibition assays.
<b>Expiration Date</b>	12 months from date of receipt.