
Product Datasheet

Anti-Influenza B Neuraminidase Antibody (orb1148910)

Description

Anti-Influenza B Neuraminidase Mouse

Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, NeA
Immunogen	19601 = Influenza B Lee/40; 19602, 19603 = Influenza B HK/73. Host Species: Mouse Antibody Class: IgG2a Stabilizers: None. Preservatives: None. Available on request.
Target	Influenza B Neuraminidase
Preservatives	PBS, pH 7.4.
Form/Appearance	Liquid
Storage	These antibodies are stable for at least one (1) year at -20 to -70°C. Store product in appropriate aliquots to avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P03474
Dilution Range	These antibodies may be used in immunoassays to detect and quantitate Influenza B neuraminidase and in neuraminidase inhibition assays.
Expiration Date	12 months from date of receipt.