

## Product Datasheet

# Anti-Canine distemper virus Antibody (orb1148922)

<b>Description</b>	Anti-Canine distemper virus Mouse Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Canine distemper virus
<b>Target</b>	Canine distemper virus
<b>Preservatives</b>	0.09% Sodium Azide. PBS, pH 7.4, 0.09% sodium azide.
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	This antibody is stable for at least one (1) year at -20C. Store in appropriate aliquots to avoid multiple freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified
<b>Dilution Range</b>	ELISA: For sandwich ELISA use #18935 as capture antibody and #18934 as detection antibody. Endusers should determine optimal antibody concentrations for their preparations of CDV.
<b>Expiration Date</b>	12 months from date of receipt.