

---

## Product Datasheet

### Anti-Bcl-xL (Phospho-Ser62) Antibody (orb1148932)

**Description**

Anti-Bcl-xL (Phospho-Ser62) Rabbit

<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IF, IHC, WB
<b>Immunogen</b>	Peptide sequence that includes phosphorylation site of serine 62 derived from human Bcl-xL and conjugated to KLH. Host Species: Rabbit Specificity: This antibody detects endogenous human, mouse, and rat Bcl-xL only when phosphorylated at serine 62. Accession no.: NP_612815.1
<b>Target</b>	Bcl-xL (Phospho-Ser62)
<b>Preservatives</b>	0.02% Sodium Azide. PBS (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1.0 mg/mL
<b>Storage</b>	This antibody is stable for at least one (1) year at -20C. Can be stored at 4C for short-term.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Purified
<b>Dilution Range</b>	Immunoblotting: use at dilution of 1:500-1:1,000. A band of ~30kDa is detected. Immunohistochemistry: use at dilution of 1:50-1:100. Immunofluorescence: use at dilution of 1:100-1:200. These are recommended working dilutions. Enduser should determine optimal dilutions for their application.
<b>Expiration Date</b>	12 months from date of receipt.