
Product Datasheet

Anti-Bcl-2 (Phospho-Ser70) Antibody (orb1148934)

Description

Anti-Bcl-2 (Phospho-Ser70) Rabbit

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Peptide sequence that includes phosphorylation site of serine 70 (R-T-S(p)-P-L) derived from human Bcl-2 and conjugated to KLH. Host Species: Rabbit Specificity: This antibody detects endogenous human Bcl-2 only when phosphorylated at serine 70. Accession no.: P10415, NP_000624.2
Target	Bcl-2 (Phospho-Ser70)
Preservatives	0.02% Sodium Azide. PBS (without Mg ² and Ca ²), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1.0 mg/mL
Storage	This antibody is stable for at least one (1) year at -20C. Can be stored at 4C for short-term.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	P10415
Entrez	151430
Dilution Range	Immunohistochemistry: use at dilution of 1:50-1:100. These are recommended working dilutions. Enduser should determine optimal dilutions for their application.
Expiration Date	12 months from date of receipt.