

Product Datasheet

Anti-Amitriptyline Antibody (orb1148966)

Description	Anti-Amitriptyline Mouse Monoclonal Antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Amitriptyline-KLH.
Target	Amitriptyline
Preservatives	PBS, pH 7.4.
Form/Appearance	Liquid
Storage	These antibodies are stable for at least one (1) year at -20C to -70C. Store product in appropriate aliquots to avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
CAS Number	50-48-6
Dilution Range	These antibodies have been qualified for use in ELISA to detect amitriptyline.
Expiration Date	12 months from date of receipt.