

Product Datasheet

Anti-Adenovirus Antibody (orb1149018)

Description	Anti-Adenovirus Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide corresponding to amino acids 255-264 of the 34 kD fiber knob protein of adenovirus 3. Host Species: Rabbit Stabilizers: None Preservatives: 0.02% sodium azide.
Target	Adenovirus
Preservatives	PBS, pH 7.4
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -20C. Avoid repeated freezing and thawing.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Dilution Range	Immunoblotting: use at 1:500-1:1,000.
Expiration Date	12 months from date of receipt.